


Thesis - Camera Calibration Instructions

Launch : `code/catkin_ws/src/camera_calibration/launch/external.launch`

Node : `code/catkin_ws/src/camera_calibration/nodes/eye_to_hand_calibration.py`



Settings

Config files can be found at

```
code/catkin_ws/src/camera_calibration/config/
```

Default Camera Configs

Eye-in-Hand Camera @ `default_eye_in_hand.json` with **large** calibration board

Top Camera @ `cam_top_default.json` with **small** board board

Front Camera @ `cam_front_default.json` with **small** board board

Example :

```
{
    "board_name": "large", // boards in boards.json
    "camera_name": "cam_front_default", // cameras in cameras.json
    "mode": "eye_to_hand", // or eye_to_hand
    "camera_topic": "/camera_front/color/image_raw", //
camera/cam_top/cam_front
    "memory_size" : 50, // max limit of camera transforms
    "load_data_directory": null, // name of file in
external_calibration_data
    "save_data_directory": "cam_top" // stored in calibration_results
}
```

Start Everything

1. Attach ChArUco in robot hand
2. Start Camera
3. Start Arm
4. Start Rviz
5. Start Calibration

```
roslaunch camera_calibration external.launch config=cam_front_default
```

Calibrate

Controls Overview

```
q = quit
c = collect transform
u = undo last transform
r = estimate pose
p = plot and publish pose estimation
s = save camera estimation
```

Collect Transforms

1. Move arm so that ChArUco is visible to the camera
2. Check quality
 - More blue dots = more stable estimation
 - The axis align with the board
3. Press the **C** key to collect transforms
4. repeat for **3** to **N** times

Publish Estimation

Press the **P** key to publish the estimation in the TF Tree

Save

Press **S** to save the pose estimation to :

```
camera_calibration/calibration_results/eye_{to/in}_hand/{camera}/{filename}
e}
```